



SURGICAL INSTRUMENT PRODUCT MANUAL



RWD Life Science Co.,Ltd

Add: Building 7A (9th, 19th-20th Floors), Building 7D (9th Floor), Shenzhen International Innovation Valley,
Dashi 1st Road, Nanshan District, Shenzhen, Guangdong, China
Web: www.rwdstco.com E-mail: rwd@rwdstco.com

RWD Life Science Inc.

Add: 10410 Corporate Drive, Sugar Land, TX 77478, USA
Tel: 1-858-900-5879 Tech support: service@rwdstco.com



Amazon

RWD Life Science

INTRODUCTION

Established in 2002, RWD Life Science has been dedicated to being a leading manufacturer in research and veterinary equipments worldwide. We specialize in anesthesia, physiological and surgical equipments, used by pre-clinical researchers and veterinarians. Our mission is to provide high quality products and integrated solutions for our customers with prompt and professional support.

For more information about our supplies, such as Laser Speckle Perfusion Imager, Anesthesia Machine, Ventilator, Stereotaxic Instrument, Syringe Pump, fiber photometry system and Surgical tools, please check our website: www.rwdstco.com

CATALOGUE

01	SURGICAL SCISSORS	02
02	BONE INSTRUMENTS	11
03	SURGICAL FORCEPS	13
04	HEMOSTATIC INSTRUMENTS	24
05	SUTURING INSTRUMENTS	29
06	SCALPELS & BLADES	33
07	RETRACTING INSTRUMENTS	35
08	CURETTING & SCRAPING INSTRUMENTS	37
09	ANIMAL IDENTIFICATION INSTRUMENTS	38
10	TISSUE MATRICES	39
11	DISINFECTION AND STERILIZATION	41
12	INSTRUMENT KIT	42

1.SURGICAL SCISSORS

Spring Scissors



01 S11001-08

- Triangular/Straight
- Sharp&Sharp/8 cm
- Edge: 8*1.55mm



02 S11002-08

- Triangular/Curved
- Sharp&Sharp/8 cm
- Edge: 8*1.55mm



03 S11035-08

- Triangular/Straight
- Sharp&Sharp/8 cm
- Edge:3.5*1.3mm



04 S11037-08

- Flat/Curved
- Sharp&Sharp/8 cm
- Edge:8*1.5mm




05 S11036-08

- Flat/Straight
- Sharp&Sharp/8 cm
- Edge:8.2*1.5mm



06 S11003-08

- Triangular/Angled
- Sharp&Sharp/8.5 cm
- Edge: 6.3*1.35mm



07 S11008-12

- Flat/Curved
- Sharp&Sharp/12cm
- Edge:14*4.5mm



08 S11007-12

- Flat/Straight
- Sharp&Sharp/12.5cm
- Edge:19*4.25mm



09 S11014-12

- Flat/Straight
- Sharp&Sharp/12.5cm
- Edge:5.5*1.9mm



10 S11019-14

- Flat/Straight
- Sharp&Sharp/14cm
- Edge:14.5*3.2mm



11 S11022-14

- Flat/Curved
- Sharp&Sharp/14cm
- Edge:14.5*3.2mm

Fine Scissors



01 S12000-09

- Sharp&Sharp/Straight
- Edge:24*4.2mm
- Length: 9cm




02 S12001-09

- Sharp&Sharp/Curved
- Edge:24*4.6mm
- m



03 S12047-10

- Sharp&Sharp/Straight
- Edge:18*4mm
- Length: 10cm



04 S12008-10

- Sharp&Sharp/Curved
- Edge:21*4.6mm
- Length: 10.5cm



05 S12007-10

- Sharp&Sharp/Straight
- Edge:21*4.2mm
- Length: 10.5cm



06 S12011-11


- Sharp&Sharp/Straight
- Edge:25*5mm
- Length: 11.5cm



07

S12012-11

- Sharp&Sharp/Straight
- Edge:25*5.4mm
- Length: 11.5cm



08

S12031-11

- Sharp&Sharp/Full Curved
- Edge:29*6mm
- Length: 11.5cm



03

S13003-11

- Blunt&Blunt/Straight
- Edge:25*4mm
- Length: 11.5cm



04

S13004-11

- Blunt&Blunt/Curved
- Edge:28*5mm
- Length: 11.5cm



09

S12011-12

- Sharp&Sharp/Straight
- Edge:31*5.8mm
- Length: 12.5cm



10

S12034-12

- Sharp&Sharp/Straight
- Edge:17*6mm
- Length: 12.5cm



05

S13020-11

- Sharp&Sharp/Straight
- Edge:20*5mm
- Length: 12cm



06

S13029-14

- Blunt&Blunt/Straight
- Edge:37*9mm
- Length: 14cm

Dissecting Scissors



01

S13001-10

- Sharp&Sharp/Straight
- Edge:25.5*6mm
- Length: 10cm



02

S13002-10

- Sharp&Sharp/Curved
- Length: 10cm




07

S13030-14

- Blunt&Blunt/Straight
- Edge:41*9.5mm
- Length: 14cm

Operating Scissors



01 S14014-10

- Sharp&Sharp/Straight
- Edge:32.7*8.4mm
- Length: 10.5cm



02 S14001-10

- Sharp&Blunt/Straight
- Edge:33*7.5mm
- Length: 11cm



03 S14016-11

- Sharp&Sharp/Curved
- Edge:35.4*7.5mm
- Length: 11.5cm



04 S14014-12

- Sharp&Sharp/Straight
- Edge:37.6*8.5mm
- Length: 12.5cm



05 S14015-12

- Blunt&Blunt/Straight
- Edge:37*7mm
- Length: 12.5cm



06 S14014-13

- Sharp&Sharp/Straight
- Edge:45*9.8mm
- Length: 13cm



07 S14016-13

- Sharp&Sharp/Curved
- Edge:40*9.2mm
- Length: 13cm



08 S14001-13

- Sharp&Blunt/Straight
- Edge:40.5*9.4mm
- Length: 13cm



09 S14001-15

- Sharp&Blunt/Straight
- Edge:50*10mm
- Length: 15.5cm



10 S14014-15

- Sharp&Sharp/Straight
- Edge:48*10mm
- Length: 15.5cm

Corneal & Cataract Scissors



01 SC11002R-11

- Corneal Scissors/11cm
- Sharp&Sharp/Curved
- Edge:11*3mm



02 SC11002T-11

- Corneal Scissors/Titanium
- Sharp&Sharp/Curved/Length: 11cm
- Edge:11*3mm




03 SC11001R-12

- Corneal Scissors/
Length: 12cm
- Sharp&Sharp/Curved
- Edge:15*3mm



04 SC11001T-12

- Corneal Scissors/Titanium/
Curved
- Sharp&Sharp/Length: 12cm
- Edge:15*3mm



03 S15006-09

- Stitch Scissors
- Sharp&Sharp/Angled
- Edge:13.5*4mm
- Length: 9.5cm



05 SC12001R-13

- Cataract Scissors/
Length: 13cm
- Sharp&Sharp/Curved
- Edge:19*1mm



06 SC12002R-13

- Cataract Scissors/
Length: 13cm
- Sharp&Sharp/Curved
- Edge:19*1mm

Other Scissors



01 S15001-09

- Ligature Scissors
- Slender Type
- Edge:20.3*4mm
- Length: 9cm



02 S15001-11

- Ligature Scissors
- Slender Type
- Edge:27*6mm
- Length: 11cm

2.BONE INSTRUMENTS

Bone Cutters

01 S22008-10



- Convex Blades
- Edge: 7.5*4.3mm
- Length: 10cm

02 S16001-10




- Straight/Length: 10.5cm
- Edge: 21*10mm

03 S22004-11




- Flat Blades/Length: 11cm
- Edge:15*6.8mm

04 S22005-12




- Flat Blades/Length: 12.5cm
- Edge:15.2*7.4mm

03 S21019-14




- Straight (SGL)
- Groove:0.5*7.5mm
- Length:14.5cm

04 S21023-14



- Curved(SGL)
- Groove: 0.7*7mm
- Length:14cm

05 S21002-15



- Medium Curved(SGL)
- Groove:3.5*12.5mm
- Length:15cm


06 S21010-15



- Light Curved(Double)
- Groove:2*19mm
- Length:15.5cm


Bone Rongeurs

01 S21022-14



- Curved(SGL)
- Groove Width:0.5mm
- Length:14cm

02 S21024-14



- Curved(SGL)
- Groove:1*6.5mm
- Length:14cm

3.SURGICAL FORCEPS

Micro Forceps

01 F11008-09



- 1x2 teeth/Straight/Length: 9.5cm
- Head length 22.5mm
- Tip 0.2*0.4mm/5mm Platform

02 F11020-11



- Straight/Length: 11cm
- Head length 30mm
- Tip 0.05*0.01mm

03 F11031-11



- Titanium/Curved/Length: 11cm
- Head length 35.5mm
- Tip 0.07*0.03mm

04 F11021-11



- 45°Angled/Length: 11cm
- Head length 30mm
- Tip 0.1*0.06mm

05 F11030-11



- Titanium/45°Angled/Length: 11cm
- Head length 32mm
- Tip 0.1*0.06mm

06 F11029-11



- Titanium/Straight/Length: 11cm
- Head length 30mm
- Tip 0.14*0.09mm

07 F11002-11



- Straight/Length: 11cm
- Head length 30mm
- Tip 0.15*0.1mm

08 F11001-11



- Straight/Length: 11cm
- Head length 30mm
- Tip 0.2*0.12mm

09 F11004-11



- Curved/Length: 11cm
- Head length 30mm
- Tip 0.2*0.12mm

10 F11003-11



- 45°Angled/Length: 11cm
- Head length 30mm
- Tip 0.3*0.15mm

11 F13002-11



- Straight/Length: 11cm
- Head length 30mm
- Tip 0.35*0.25mm

12 F11023-11



- Straight/Length: 11cm
- Head length 29mm
- Tip 0.1*0.06mm

13 F13004-11



- Curved/Length: 11.5cm
- Head length 30mm
- Tip 0.3*0.2mm

14 F11022-11



- Curved/Length: 11.5cm
- Head length 32mm
- Tip 0.1*0.06mm

15 F11027-13



- Straight/Length: 13cm
- Head length 37.5mm
- Tip 0.1*0.06mm

16 F11028-13



- Curved/Length: 13cm
- Head length 38mm
- Tip 0.1*0.06mm

17 F11001-13



- Straight/Length: 13cm
- Head length 33mm
- Tip 0.15*0.1mm

18 F11003-13



- 45°Angled/Length: 13cm
- Head length 30mm
- Tip 0.15*0.1mm

Dissecting Forceps

01 F14012-10



- Light Curved
- Tip 13*0.5mm
- Length: 10cm

02 F12005-10



- Straight/10cm
- Tip 14*1mm
- Length: 10cm

03 F14013-10



- Large Curved
- Tips Width 0.5 mm
- Length: 10cm

04 F12007-10



- Large Curved
- Tip 10*1mm
- Length: 10cm

05 F12006-10



- Light Curved
- Tip 13.5*1mm
- Length: 10cm

06 F12029-09



- Straight/Length: 9.5cm
- Tip 13*0.8mm

07 F12012-10



- Straight/Length: 10cm
- Tip 0.4*0.2mm

08 F12013-10



- 45°Angled/Length: 10cm
- Tip 0.5*0.2mm

13 F12003-12



- Length: 12cm
- Tips width: 1.5mm

14 F12011-13



- Curved/Length: 13cm
- Tip width: 1.2mm

09 F12010-10



- Straight/Length: 10.5cm
- Tip 15*1.8mm

10 F12019-10



- Straight/Length: 10.5cm
- Tip 15*2mm

15 F12010-13



- Straight/Length: 13cm
- Tip 17*1.3mm

16 F15007-15



- Angled/Length:15cm
- Tip width 1.2mm

11 F12011-10



- Curved/Length: 10.5cm
- Tip 15*1.9mm

12 F12002-12



- Length: 12cm
- Tip 11.5*0.6mm

17 F12016-15



- Straight/Length: 15.5cm
- Tip 19*1.4mm

Tissue Forceps(Toothed)

01 F13005-08




- 1x2 Teeth/Length: 8.5cm
- Tip width: 0.7mm

02 F13029-10



- 1x2 Teeth/Straight
- Tip width 1.3mm
- Length: 10.5cm

03 F13030-10



- 1x2 Teeth/Light Curved
- Tip width 1.3mm
- Length: 10.5cm

04 F13019-12



- 1x2 Teeth/Length: 12cm
- Tip width 1.8mm

05 F13023-14



- 1x2 Teeth/Straight
- Tip width 1.8mm
- Length: 14cm

Corneal Forceps

01 FC12002R-08



- Arc/Length:8.5cm
- Tip 28*0.3mm/6mm platform

02 FC12002T-08



- Titanium/Arc
- Tip 28*0.3mm/6mm platform
- Length: 8.5m

03 FC13002R-08



- 1×2 Teeth/Straight
- Tip 0.3*0.5mm/6mm platform
- Length: 8.5m

04 FC13002T-08



- Titanium/1×2 Teeth/Straight
- Tip 0.3*0.5mm/6mm platform
- Length: 8.5m

05 FC12001R-11



- Straight/Length: 10.7cm
- Tip 28*0.3mm/6mm platform

06 FC12001T-11



- Titanium/Straight
- Tip 28*0.3mm/6mm platform
- Length: 10.7cm

07 FC13001R-11

- 1×2 Teeth/Straight
- Tip 0.3*0.5mm/6mm platform
- Head length: 28cm

08 FC13001T-11

- Titanium/1×2 Teeth/Straight
- Tip 0.3*0.5mm/6mm platform
- Head length: 28cm

Tying Forceps

01 FC11002R-08

- Straight/Length: 8.5cm
- Head length 26.5mm/
Tip 0.3*0.5mm
- 6mm Tying platform

02 FC11002T-08

- Titanium/Straight/
Length: 8.5cm
- Head length 26.5mm/
Tip 0.3*0.5mm
- 6mm Tying platform

03 FC11001R-11

- Straight/Length: 11cm
- Head length 23.5mm/
Tip 0.3*0.5mm
- 6mm Tying platform

04 FC11001T-11

- Titanium/Straight/
Length: 11cm
- Head length 23.5mm/
Tip 0.3*0.5mm
- 6mm Tying platform

Capsulorhexis Forceps

01 FC14001R-11

- Arc/Round handle
- Head length 14mm/
Tip 0.6*0.6mm
- Length: 11cm

02 FC14001T-11

- Titanium/Arc/
Round handle
- Head length 14mm/
Tip 0.6*0.6mm
- Length: 11cm


Other Forceps

01 F18011-09


- Suture Forceps
- Teeth/Straight
- Tip 0.4*0.9mm/
5mm Platform
- Length: 9cm

02 F18017-09

- Ring Forceps
- OD×ID 3.2×2.2mm
- Length: 9cm

03 F18015-10

- Cover Glass Forceps
- Straight /Length: 10.5cm

04 F18001-10

- Cover Glass Forceps
- Angled/Length: 10.5cm
- Tip Size 12*4.7mm

05 F16006-11



- Towel Forceps
- Length: 11cm

06 F18008-15



- Tissue Forceps
- 4x5 Teeth/Length:15cm
- Tip 39*5mm

4.HEMOSTATIC INSTRUMENTS

Hemostatic Forceps

01 F22002-10



- Straight /Length:10.5cm
- Tip: 20*0.8mm

02 F22003-10



- Curved/Length:10cm
- Tip: 20.8*1mm

03 F21001-12



- Straight /Length:12.5cm
- Tip: 23.5*1mm

04 F21002-12



- Curved/Length:12.5cm
- Tip:24*1mm

05 F22010-12




- Straight /Length:12.5cm
- Tip width 1.2mm

06 F21003-12



- 1x2 Teeth
- Tip: 18*1.5mm
- Length:12.5cm

07 F21011-14



- Straight /Length:14.5cm
- Tip: 36*2mm

08 F21020-16




- Curved/Length:16cm
- Tip: 35*1.5mm

09 F21015-16



- Straight /Length:16cm
- Tip: 33*1.8mm

10 F21019-16



- Straight /Length:16cm
- Tip: 37*1mm

01 R33002-38



- Curved/Length:38mm
- Tip: 11*3.8mm

02 R33001-48



- Straight /Length:48mm
- Tip: 14*4mm

03 R33002-48



- Curved/Length:48mm
- Tip:14*4mm

04 R33003-02



- Curved/Length:58mm
- Tip:18*4mm

05 R33003-01



- Straight /Length:58mm
- Tip: 18*4mm

06 R31010-26



- 45°Angled/Length:26.5mm
- Tip: 15*2.5mm

07 R31009-26



- 30°Angled/Length:28mm
- Tip: 14*2.5mm

08 R31008-26



- Straight /Length:28mm
- Tip: 11*2mm

09 R34001-14



- Length: 14 cm
- Tip: 8.4*5mm
- For R31005- and R31006- Clamps

10 R31006-06



- Curved/Length: 16mm
- Tip: 4*0.75mm

15 R31005-08



- Straight /Length: 18mm
- Tip: 8*2mm

11 R31006-04



- Curved/Length: 17mm
- Tip: 6*1mm

12 R31005-06

- Straight /Length: 19mm
- Tip: 6*1mm

13 R31005-04



- Straight /Length: 19mm
- Tip: 16mm

14 R31005-10



- Straight /Length: 20mm
- Tip: 10*2mm

5.SUTURE INSTRUMENTS

Needle Holders

01 F31047-12



- Straight /Length:12cm
- Tip: 8.5*1.5mm

02 F31017-12




- Straight /lock & w/o saws
- Length: 12 cm

03 FC21001R-12



- Curved/Length:12.5cm
- Tip width 0.35mm

04 F31024-13



- Straight /w/o Serrations
- Length: 13 cm

05 F31025-13



- Straight /Length:13cm
- Tip: 17.5*1.75mm

06 F31047-14



- Straight /Length:14cm
- Tip: 17.3*1.7mm

07 F31029-15



- Straight /Length:15cm

08 F31031-15



- Straight /Length:15.5cm
- Tip: 21*2.5mm

Suture Needle



F33304-10	○1/2/4×10mm	50/Box
F33306-14	○1/2/6×14mm	50/Box
F33307-17	○1/2/7×17mm	50/Box
F33309-24	○1/2 9×24mm	50/Box



F33104-10	$\triangle 1/2/4\times 10\text{mm}$	50/Box
F33105-12	$\triangle 1/2/5\times 12\text{mm}$	50/Box
F33106-14	$\triangle 1/2/6\times 14\text{mm}$	50/Box
F33107-17	$\triangle 1/2/7\times 17\text{mm}$	50/Box


Suture



F35403-30	$\bigcirc 3/8$ Needle/ $5\times 14\text{mm}$	Size 3-0/90cm	50/Box
F35401-50	$\bigcirc 3/8$ Needle/ $3\times 10\text{mm}$	Size 5-0/90cm	50/Box
F35204-70	$\triangle 3/8$ Needle/ $2.5\times 5\text{mm}$	Size 7-0/30cm	50/Box
F35207-30	$\triangle 3/8$ Needle/ $5\times 14\text{mm}$	Size 7-0/90cm	50/Box



F34004-20	Size 2-0/60cm	15Pcs/Pkg	40 Pkg/Box
F34004-30	Size 3-0/60cm	15Pcs/Pkg	40 Pkg/Box
F34004-40	Size 4-0/60cm	15Pcs/Pkg	40 Pkg/Box
F34004-50	Size 5-0/60cm	15Pcs/Pkg	40 Pkg/Box
F34004-60	Size 6-0/60cm	15Pcs/Pkg	40 Pkg/Box
F34004-00	Size 0#/60cm	15Pcs/Pkg	40 Pkg/Box




18020-20	Size 2/0 / 0.32mm	Length: 22.5cm
18020-40	Size 4/0 / 0.17mm	Length: 22.5cm
18020-60	Size 6/0 / 0.08mm	Length: 22.5cm

6.SCALPELS & BLADES

Scalpel

01 S33007-12



- Cutting Edge: 8mm
- Width: 1.5mm
- Length: 12 cm

02 S33006-13



- Length: 13 cm
- Tip: 25*3mm

Scalpel Blades

01 S31010-01



- 10# Scalpel Blades
- Box of 100pcs

02 S31011-01



- 11# Scalpel Blades
- Box of 100pcs

03 S31012-01



- 12# Scalpel Blades
- Box of 100pcs

04 S31015-01



- 15# Scalpel Blades
- Box of 100pcs

05 S31021-02



- 21# Scalpel Blades
- Box of 100pcs

06 S31022-02



- 22# Scalpel Blades
- Box of 100pcs

07 S31023-02



- 23# Scalpel Blades
- Box of 100pcs

01 S33009-06



- Scalpel Blade Remover
- Length: 6.5cm

Scalpel Handles

01 S32003-12



- Length: 12.5cm


02 S32004-13



- Length: 14cm

7.RETRACTING INSTRUMENTS

Retractors

A U-shaped metal retractor with two loops at the top.


01 R22029-03

- 11mm spread
- Length: 3cm

A U-shaped metal retractor with two loops at the top, wider than R22029-03.

02 R22029-04

- 20mm spread
- Length: 4cm

A retractor with two long, curved arms and multiple teeth.


03 R22005-45

- 3x3 Teeth/Blunt
- 36mm Spread
- Length: 5cm

A retractor with two long, curved arms and multiple teeth, similar to R22005-45 but with a different spread.

04 R22009-01

- 4x4 Teeth/Blunt
- 65mm Spread
- Length: 7.5cm

A retractor with two long, curved arms and multiple teeth, similar to R22005-45 but with a different spread.

05 R22013-10

- 2x3 Teeth/Blunt
- Length: 10cm

Hooks

A long, thin metal hook with a curved tip.


01 R21015-16

- 45 ° Curved
- 1 Ball Tooth/Depth 10mm
- Length: 16cm

A long, thin metal hook with a curved tip.

02 R42002-12


- Tip Dia 3mm
- Length: 12cm

A long, thin metal hook with a curved tip.

03 R21013-01

- 1 Tooth/Blunt
- Length: 15.5cm

Micro Spatula

A long, thin metal spatula with a curved tip.

01 R51003-11

- Curved/Length:11cm
- Tip:2*0.3cm

8.SCRAPING / CURETTING INSTRUMENTS

Spoons



01 R52001-12

- Dia. 8 mm x Dep. 2 mm
- Length: 12.5cm

Chisels



01 S23002-16

- Curved/Length:16cm
- Tip:12*4cm

9.ANIMAL IDENTIFICATION INSTRUMENTS

Finger Loop Ear Punch

01 C22001-10



- Dia. 1 mm
- Length: 10.5cm

02 C22002-10



- Dia. 2 mm
- Length: 10.5cm

Ear Tag applicator

01 23001-5



- Dia. 2.55 mm
- Length: 14.5cm

Ear Tag

01 C23001-14



- 9.2*2.5*0.4mm

Brain Matrix

10.TISSUE MATRIX


68707



68709



68711



68712



68713



Model	Applicable Animals	Direction
68707	Mouse/20-75g	Coronal/1mm
68708	Mouse/20-75g	Sagital/1mm
68709	Rat/75-300g	Coronal/1mm
68710	Rat/75-300g	Sagital/1mm
68711	Rat/300-600g	Coronal/1mm
68712	Rat/300-600g	Sagital/1mm
68713	Mouse/20-75g	Coronal/0.5mm

Heart Matrix

68719



68720



68721



68722




Model	Applicable Animals	Direction
68719	Mouse/Neonatal Rat/20-75g	Coronal/1mm
68720	Rat/175-400g	Coronal/1mm
68721	Mouse/Neonatal Rat/20-75g	Coronal/0.5mm
68722	Rat/175-400g	Coronal/0.5mm

11.DISINFECTION AND STERILIZATION


Sterilization Trays

01 C12006-01



• 270*160*22mm

02 C12006-02



• 215*135*22mm

Sterilizer

01 RS1500



• Weight: 1500g

02 RS3002



- Volume: φ4x14cm
- Power: 220V

12.INSTRUMENT KIT

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
General Surgery Instrument Kit for Rodents	SP0001-G	VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		IRIS-Fine Scissors-S/S Str/11.5cm	S12011-11	1
		IRIS-Fine Scissors-S/S Cvd/11.5cm	S12012-11	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/9cm	S12001-09	1
		Micro Forceps-Str,Tip 0.2*0.12mm, 11cm	F11001-11	1
		Micro Forceps-Cvd,Tip 0.2*0.12mm, 11cm	F11004-11	1
		Dressing Forceps-Str, Tip width 1.8mm,10.5cm	F12010-10	1
		Dressing Forceps-Cvd, Tip width 1.9mm,10.5cm	F12011-10	1
		HALSTED Artery Forceps-Str/1mm/12.5cm	F21001-12	2
		HALSTED Artery Forceps-Cvd/1mm/12.5cm	F21002-12	2
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	2
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	2
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		MAYO-STILLE Dissecting & Operating Scissors-B/B Str/14cm	S13029-14	1
		SPENCER Ligature Scissors (Slender Type)-11cm	S15001-11	1
		BOHLER Rongeurs (Double)-Light Cvd/2mm/15.5cm	S21010-15	1
		Bone Cutters with Flat Blades (SGL)/11cm	S22004-11	1
		Operating Scissors (Round Type)-S/S Str/12.5cm	S14014-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		Scalpel Blade Remover-6.5cm	S33009-06	1
		ADSON 1x2 Teeth Tissue Forceps, Tip width 1.8mm, 12cm	F13019-12	1
		MANNERFELT Hooks, 1 Tooth, Blunt, 15.5cm	R21013-01	1
		ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm	R22009-01	1
		3x3 Teeth Retractors-Blunt, 30mm Spread, 5cm	R22005-45	1
		Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	R51003-11	1
		NA Nylon Sutures w/Needle-△3/8/7-0 (50/Box)	F35204-70	20
		NA Terylene Sutures w/Needle-○3/8/ 5-0 (50/Box)	F35401-50	20
		SS Micro Clamps-Cvd/L*W 4*0.75mm/16mm	R31006-06	5
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	5
		Clip Applicator for R31005- and R31006-Clamps/14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
General Ophthalmic Surgery Instrument Kit for Rodents	SP0002-G	VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		VANNAS Spring Scissors (Triangular)-S/S Angled/8.5cm	S11003-08	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/10.5cm	S12008-10	1
		O'BRIEN Stitch Scissors-S/S Angled/9.5cm	S15006-09	1
		Dressing Forceps w/o Serrations-45°Cvd, S/S, Tip 0.5*0.2mm, 10cm	F12013-10	1
		ADSON 1x2 Teeth Tissue Forceps, Tip width 1.8mm, 12cm	F13019-12	1
		IRIS Dissecting Forceps-Str, Tip width 1mm, 10cm	F12005-10	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, 10cm	F12006-10	1
		McPHERSON Forceps w/5mm Platform-Str, Tip 0.4*0.9mm, 9cm	F18011-09	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		MANNERFELT Hooks, 1 Tooth, Blunt, 15.5cm	R21013-01	1
		Colibri Eye Specula-20mm spread, 4cm	R22029-04	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Microsurgery Instrument Kit for Rat, Cat, Rabbit etc.	SP0003-R	IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/11.5cm	S12012-11	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		Micro Forceps-Str, head length 30mm, Tip 0.2*0.12mm, 11cm	F11001-11	1
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	5
		Clip Applicator for R31005- and R31006-Clamps/14cm	R34001-14	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		HALSTED Mosquito Forceps-Str, 1mm, 12.5cm	F22001-12	1
		HALSTED Mosquito Forceps-Cvd,1mm, 12.5cm	F22002-12	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Microsurgery Instrument Kit for Mice	SP0003-M	IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/9cm	S12001-09	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		Micro Forceps-Str, Tip 0.2*0.12mm, 11cm	F11001-11	1
		SS Micro Clamps-Str/L*W 4*0.75mm/16mm	R31005-04	5
		Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm	R34001-14	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	1
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Cardio-cerebral Perfusion Surgery Instrument Kit for Mice	SP0008-M	Operating Scissors (Round Type)-S/S Str/10.5cm	S14014-10	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		IRIS Dissecting Forceps-Large Cvd, Tip width 1mm, 10cm	F12007-10	1
		Dressing Forceps-Str, Tip width 1.8mm,10.5cm	F12010-09	2
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	2
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	2
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		NA Terylene Sutures w/Needle-△3/8 /7-0 (50/Box)	F35204-70	20
		MANNERFELT Hooks, 1 Tooth, Blunt, 15.5cm	R21013-01	1
		3x3 Teeth Retractors-Blunt, 36mm Spread, 5cm	R22005-45	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		SS Micro Clamps-Str/L*W 4*0.75mm/16mm	R31005-04	5
		Clip Applicator for R31005- and R31006-Clamps-14cm	R34001-14	1
		Spinal Cord Hook - Tip Dia. 3mm/12cm	R42002-12	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Brain Microinjection Surgery Instrument Kit for Rat	SP0005-R	Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		IRIS-Fine Scissors-S/S Str, 10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd, 10.5cm	S12008-10	1
		Dressing Forceps-Str, Tip width 1.8mm,10.5cm	F12010-10	2
		Dressing Forceps-Cvd, Tip width 1.9mm,10.5cm	F12011-10	2
		HALSTED Mosquito Forceps-Str, 1mm/12.5cm	F21001-12	2
		HALSTED Mosquito Forceps-Cvd, 1mm/12.5cm	F21002-12	2
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		Operating Scissors (Round Type)-S/S Str/12.5cm	S14014-12	1
		NA Nylon Sutures w/Needle-△3/8/7-0 (50/Box)	F35204-70	20
		ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm	R22009-01	1
		Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	R51003-11	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Univeral Dissecting Surgery Instrument Kit for Mice	SP0007-M	IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/9cm	S12001-09	1
		STEVENS Fine Dissecting Scissors (Flat Type)-B/B Str/11.5cm	S13003-11	1
		STEVENS Fine Dissecting Scissors (Flat Type)-B/B Cvd/11.5cm	S13004-11	1
		SPENCER Ligature Scissors (Slender Type)-11cm	S15001-11	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		15# Scalpel Blades (Box of 100pcs)	S31015-01	1
		IRIS Dissecting Forceps-Str, Tip width 1mm, 10cm	F12005-10	2
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm,10cm	F12006-10	2
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	2
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	2
		IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm	F13029-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35305-50	20
		3x3 Teeth Retractors-Blunt, 36mm Spread, 5cm	R22005-45	1
		SS Micro Clamps-Str/L*W 4*0.75mm/16mm	R31005-04	5
		Spinal Cord Hook - Tip Dia. 3mm/12cm	R42002-12	1
		Clip Applicator for R31005- and R31006-Clamps-8.4*5mm/14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Univeral Dissecting Surgery Instrument Kit for Rat	SP0007-R	IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/10.5cm	S12008-10	1
		IRIS-Fine Scissors-S/S Cvd/11.5cm	S12012-11	1
		MAYO-STILLE Dissecting & Operating Scissors-B/B Str/14cm	S13029-14	2
		MAYO-STILLE Dissecting & Operating Scissors-B/B Cvd/14cm	S13030-14	2
		HALSTED Artery Forceps-Str/1mm/12.5cm	F21001-12	2
		HALSTED Artery Forceps-Cvd/1mm/12.5cm	F21002-12	2
		Dressing Forceps-Str, Tip width 1.8mm, 10.5cm	F12010-10	1
		Dressing Forceps-Cvd, Tip width 1.9mm, 10.5cm	F12011-10	1
		IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm	F13029-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str, 1.5mm/12cm	F31047-12	1
		SPENCER Ligature Scissors (Slender Type)/11cm	S15001-11	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		15# Scalpel Blades (Box of 100pcs)	S31015-01	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35305-50	20
		ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm	R22009-01	1
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	5
		Spinal Cord Hook - Tip Dia. 3mm/12cm	R42002-12	1
		Clip Applicator for R31005- and R31006-Clamps-14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Brain Microinjection Surgery Instrument Kit for Mice	SP0005-M	Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd/9cm	S12001-09	1
		Operating Scissors (Round Type)-S/S Str/10.5cm	S14014-10	1
		IRIS Dissecting Forceps-Str, Tip width 1mm, 10cm	F12005-10	2
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, 10cm	F12006-10	2
		HARTMAN Mosquito Forceps-Str/20*0.8mm/10.5cm	F22002-10	2
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		NA Nylon Sutures w/Needle-△3/8/2.5×5/30cm/7-0 (50/Box)	F35204-70	20
		HARTMAN Mosquito Forceps-Cvd/Tip width 1mm/10cm	F22003-10	2
		Dual-end Micro Spatula w/tip W. 2mm x Thk. 0.3mm - Cvd/11cm	R51003-11	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Cardio-cerebral Perfusion Surgery Instrument Kit for Rat	SP0008-R	Operating Scissors (Round Type)-S/S Str/12.5cm	S14014-12	1
		IRIS-Fine Scissors-S/S Str/Head length 21mm/10.5cm	S12007-10	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		IRIS Dissecting Forceps-Large Cvd, Tip width 1mm,10cm	F12007-10	1
		IRIS Dissecting Forceps-Str, Tip width 1mm,10cm	F12005-10	2
		HALSTED Mosquito Forceps-Str, 1mm/12.5cm	F21001-12	2
		HALSTED Mosquito Forceps-Cvd, 1mm/12.5cm	F21002-12	2
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35305-50	20
		NA Terylene Sutures w/Needle-△3/8 /7-0 (50/Box)	F35204-70	20
		MANNERFELT Hooks, 1 Tooth, Blunt, 15.5cm	R21013-01	1
		3x3 Teeth Retractors-Blunt, 36mm Spread, 5cm	R22005-45	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		SS Micro Clamps-Str/6*1mm/16mm	R31005-06	2
		Spinal Cord Hook - Tip Dia. 3mm/12cm	R42002-12	1
		Clip Applicator for R31005- and R31006-Clamps-14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
General Surgery Instrument Kit for Mouse	SP0009-M	Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm	F11020-11	1
		Micro Forceps-45°Angled, head length 30mm, Tip 0.1*0.06mm, 11cm	F11021-11	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		Colibri Eye Specula-20mm spread, 4cm	R22029-04	1
		IRIS Dissecting Forceps-Str, Tip width 1mm,10cm	F12005-10	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm	F12006-10	1
		IRIS Dissecting Forceps-Large Cvd, Tip width 1mm, teeth length 10mm, 10cm	F12007-10	1
		IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm	F13029-10	1
		IRIS-Fine Scissors-S/S Str, 10.5cm	S12007-10	1
		IRIS-Fine Scissors-S/S Cvd,9cm	S12001-09	1
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	1
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
General Surgery Instrument Kit for Rat	SP0009-R	Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm	F11020-11	1
		Micro Forceps-45°Angled, head length 30mm, Tip 0.1*0.06mm, 11cm	F11021-11	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		11# Scalpel Blades (Box of 100pcs)	S31011-01	1
		NA Terylene Sutures w/Needle-○3/8 /5-0 (50/Box)	F35401-50	20
		NA Terylene Sutures w/Needle-△3/8 /7-0 (50/Box)	F35204-70	20
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.5mm/12cm	F31047-12	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Cvd/8cm	S11002-08	1
		ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm	R22009-01	1
		IRIS Dissecting Forceps-Str, Tip width 1mm,10cm	F12005-10	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm	F12006-10	1
		IRIS Dissecting Forceps-Large Cvd, Tip width 0.9mm, teeth length 10mm, 10cm	F12007-10	1
		IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm	F13029-10	1
		IRIS-Fine Scissors-S/S Str/ 11.5cm	S12011-11	1
		IRIS-Fine Scissors-S/S Cvd/11.5cm	S12012-11	1
		HALSTED Artery Forceps-Str/1mm/12.5cm	F21001-12	1
		HALSTED Artery Forceps-Cvd/1mm/12.5cm	F21002-12	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Basic Surgical Instrument Kit for Rat & Mice	SP0010-B	VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS Dissecting Forceps-Str, Tip width 1mm, 10cm	F12005-10	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm	F12006-10	1
		Micro Forceps-Str, head length 36mm, Tip 0.05*0.01mm, 11cm	F11020-11	1
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	1
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str, 1.5mm/12cm	F31047-12	1
		ZIEGLER Knives-8mm Cutting Edge/1.5mm Width/12cm	S33007-12	1
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	1
		SS Micro Clamps-Cvd/L*W 4*0.75mm/16mm	R31006-06	1
		Clip Applicator for R31005- and R31006-Clamps/14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Lite Surgical Instrument Kit	SP0010-L	VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		IRIS Dissecting Forceps-Str, Tip width 1mm, 10cm	F12005-10	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm, teeth length 13.5mm, 10cm	F12006-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str, 1.5mm/12cm	F31047-12	1
		Sterilization Trays-215*135*22mm	C12006-02	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
MCAO Surgical Instrument Kit for Rat & Mice	SP0011-B	Micro Forceps-Str,Tip 0.1*0.06mm, 13cm	F11027-13	1
		Scalpel Handles with Ruler 3# Solid-12.5cm	S32003-12	1
		IRIS Dissecting Forceps-Light Cvd, Tip width 1mm,10cm	F12006-10	1
		VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		IRIS 1x2 Teeth Tissue Forceps-Light Cvd, Tip width 1.3mm, 10.5cm	F13030-10	1
		HARTMAN Mosquito Forceps-Str/0.8mm/10.5cm	F22002-10	1
		IRIS-Fine Scissors-S/S Cvd/10.5cm	S12008-10	1
		NA Terylene Sutures w/Needle-△3/8 /7-0 (50/Box)	F35204-70	20
		Operating Scissors (Round Type)-S/S Str/15.5cm	S14014-15	1
		MANNERFELT Hooks, 1 Tooth, Blunt, 15.5cm	R21013-01	4
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str,1.7mm/14cm	F31047-14	1
		15# Scalpel Blades (Box of 100pcs)	S31015-01	1
		IRIS Dissecting Forceps-Str, Tip width 1mm,10cm	F12005-10	1
		IRIS 1x2 Teeth Tissue Forceps-Str, Tip width 1.3mm, 10.5cm	F13029-10	1
		Micro Forceps-Cvd,Tip 0.1*0.06mm, 13cm	F11028-13	1
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	5
		Micro Forceps-Cvd,Tip 0.2*0.12mm, 11cm	F11004-11	1
		Micro Forceps-Str,Tip 0.2*0.12mm, 11cm	F11001-11	1
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	1
		SS Micro Clamps-Cvd/L*W 6*1mm/17mm	R31006-04	1
		Clip Applicator for R31005- and R31006-14cm	R34001-14	1
		Instrument Storage Portfolio	SP0001-P	1

Kit Name	Kit Model	Instrument Description	Kit Model	Quantity
Surgical Instruments Kit DEMO(Research Version)	SP0001-DR	VANNAS Spring Scissors (Triangular)-S/S Str/8cm	S11001-08	1
		VANNAS Spring Scissors (Triangular)-S/S Angled/8.5cm	S11003-08	1
		IRIS-Fine Scissors-S/S Str/10.5cm	S12007-10	1
		Friedman-Pearson Rongeurs (SGL)-Str/0.5mm/14.5cm#	S21019-14	1
		Micro Forceps-Str, Tip 0.05*0.01mm, 11cm	F11020-11	1
		Micro Forceps (Titanium)-Cvd,Tip 0.07*0.03mm, 11cm	F11031-11	1
		Micro Forceps-Str, Tip 0.2*0.12mm, 11cm	F11001-11	1
		ADSON 1x2 Teeth Tissue Forceps, Tip width 1.8mm, 12cm	F13019-12	1
		Cover Glass Forceps-Angled/Tip Size /10.5cm	F18001-10	1
		HARTMAN Mosquito Forceps-Cvd/1mm/10cm	F22003-10	1
		OLSEN-HEGAR Needle Holders with Scissors-Str, 1.5mm/12cm	F31047-12	1
		ALM 4x4 Teeth Retractors-Blunt, 65mm Spread, 7.5cm	R22009-01	1
		SS Micro Clamps-Str/L*W 4*0.75mm/16mm	R31005-04	1
		SS Micro Clamps-Str/L*W 6*1mm/19mm	R31005-06	1
		Clip Applicator for R31005- and R31006-Clamps/14cm	R34001-14	1
		STEVENS Fine Dissecting Scissors (Flat Type)-B/B Str/11.5cm	S13003-11	1
		Corneal Forceps(Titanium)-Arc, Head length 28mm, 0.3mm Tips, 6mm tying platform with φ0.2 hole, Three-hole handle, 8.5cm	FC12002T-08	1
		Instrument Storage Portfolio	SP0001-P	1

Copyright ©RWD Life Science Co., LTD. 2024. All rights reserved.

No part of this document maybe reproduced or transmitted in any form or by any means without prior written consent of RWD Life Science Co., LTD.

Trademark Notice

RWD, **RWD** | 瑞沃德, and **RWD** are trademarks or registered trademarks of RWD Life Science Co., LTD.

Other trademarks, product, service and company names mentioned are the property of their respective owners.

General Disclaimer

The information in this document may contain predictive statements including, without limitation, statements regarding the future financial and operating results, future product portfolio, new technology, etc. There are a number of factors that could cause actual results and developments to differ materially from those expressed or implied in the predictive statements. Therefore, such information is provided for reference purpose only and constitutes neither an offer nor an acceptance.

RWD may change the information at any time without notice.

